



# New Accessories for the V-600 Series UV-Vis/NIR Spectrometers

## SAH-769 One Drop Accessory



### Precise analysis of micro-volume samples for Life Science Studies

The SAH-769 One Drop accessory is a dedicated accessory for the V-600 Series to measure micro-volume samples of protein and nucleic acid. The sampling procedure is quite easy, by simply placing a droplet of a sample on the disk cell window, providing rapid measurements within 20 seconds per sample. The minimum sample volume is 5  $\mu\text{L}$  for the standard 1 mm pathlength disk cell, while 0.6  $\mu\text{L}$  can be used for the optional 0.2 mm pathlength disk cell. The shorter optical pathlength allows measurements of high concentration samples without dilution. A precise optical pathlength is maintained, assuring excellent reproducibility of measurements.



**Dedicated disk cell for SAH-769**

Optical pathlength : 1 mm (left), 0.2 mm (right)  
Sample volume : 5  $\mu\text{L}$  (left), 0.6  $\mu\text{L}$  (right)

### Quick, but highly accurate measurement by a research-grade spectrophotometer Only 20 seconds per sample.



**Step 1**

Drop a sample on the cell.




**Step 2**

Set the cover glass.



**Step 3**


Push the start button.




**Step 4**

Wipe the sample off.

### Quantitative analysis of calf thymus DNA using a cell with a 0.2 mm pathlength Only 0.6 $\mu\text{L}$ sample volume is required.



**Spectra of calf thymus DNA solutions**



**Calibration curve of calf thymus DNA**

The figure illustrates good linearity up to approx. 3000 ng/ $\mu\text{L}$ . This is equivalent to a linearity up to 35 Abs by using a 10 mm pathlength cell.

**medson**  
APARATURA ANALITYCZNA

ul. Średzka 15, 62-021 Paczkowo  
+48 61 815 74 51, NIP: PL 7770043517  
medson@medson.pl, www.medson.pl

**JASCO**




# New Accessories for the V-630Bio UV-Vis Spectrophotometer

## Dedicated Functionality for Life Science Applications


### Capillary Cell and Adapter for Melting Analysis

**Capillary cell with minimum sample volume of 3µL and 0.5 mm pathlength**




The capillary cell adapter provides temperature control of an extremely small amount of liquid samples when placed into a thermostated cell holder. The capillary cell is made of quartz with an approximately 0.5 mm optical pathlength and a minimum sample volume of 3 µL. The capillary adapter can be used with the Peltier thermostated cell holders and the 6-position cell block of the PAC-743/743R automatic 6/8-position Peltier cell changer. An optional temperature sensor for the cell block is recommended for monitoring the temperature.

#### ■ DNA melting analysis with capillary cell and 10 mm cell



**Melting curves of DNA sample measured by a capillary cell and a 10 mm cell**

The figure illustrates the DNA melting analysis of a DNA sample by using a capillary cell and a 10 mm pathlength cell with the ETCS-761 Peltier thermostated single cell holder. The melting temperature calculated from the data measured with the capillary cell was 63.91°C, while 63.81°C was obtained with the 10 mm pathlength cell. The results of both measurements agree quite well.



#### V-630Bio

Stand-alone UV-Vis Spectrophotometer

**Built-in bio application programs include:**

- Protein/nucleic acid measurement
- Temperature ramping/DNA melting analysis
- Kinetic measurements and analysis software



#### PAC-743/743R

6/8-position automatic Peltier thermostated cell changers and 6-position cell block



#### ETCS-761

Peltier thermostated single cell holder

**medson**

APARATURA ANALITYCZNA

ul. Średzka 15, 62-021 Paczkowo  
+48 61 815 74 51, NIP: PL 7770043517  
medson@medson.pl, www.medson.pl

**JASCO**



# New Accessories for the V-600 Series UV-Vis/NIR Spectrometers


## PAC-743/743R Automatic 6/8-Position Peltier Cell Changer




### PAC-743/743R Automatic 6/8-Position Peltier Cell Changer

The PAC-743/PAC-743R Automatic 6/8-Position Peltier Cell Changers allow measurements of the transmittance/absorbance of multiple sample cells by using dedicated cell blocks. With a temperature range from 0 to 100°C, four dedicated cell blocks are available, an 8-position 1 mm micro-cell block, an 8-position 10 mm micro-cell block, an 8-position 5 mm cell block and a 6-position 10 mm cell block.

### Parallel melting analysis of 10 $\mu$ L samples by using PAC-743 and 8-position 1 mm micro-cell




### Options for PAC-743 and PAC-743R




- ❶ 6-position cell block
- ❷ Capillary cell and adapter (for V-630/V-630Bio only)  
light path length: approx. 0.5 mm, minimum sample volume: 3  $\mu$ L
- ❸ Rectangular cell  
light path length: 10 mm, light path width: 10 mm
- ❹ Rectangular micro-cell  
light path length: 10 mm, light path width: 2 mm
- ❺ 8-position cell block
- ❻ Rectangular cell  
light path length: 5 mm, light path width: 5 mm
- ❼ 8-position 10 mm micro-cell block
- ❽ 8-position 10 mm micro-cell  
light path length: 10 mm, cell volume: 100  $\mu$ L

8-position 1 mm  
micro-cell block



### 8-position 1 mm micro-cell

The 8-position 1 mm micro-cell block allows measurements of a maximum of 8 samples with a 10  $\mu$ L sample volume using 1 mm optical pathlength cells. The silicon cap and cap fixture are available to prevent volatilization of samples during DNA melting measurements at high temperatures. An optional cell temperature sensor is recommended for accurate temperature readings of samples.



**medson**

APARATURA ANALITYCZNA

ul. Średzka 15, 62-021 Paczkowo  
☎ +48 61 815 74 51, NIP: PL 7770043517  
✉ medson@medson.pl, www.medson.pl

**JASCO**